

Virtual
Seminar
Shawn Lewenza

Athabasca University
University of Calgary
Luminous BioSolutions

11 am MT, March 5, 2026

Biosensors are a biological detection method for simple, rapid, specific, and quantitative naphthenic acids monitoring



ABSTRACT

Bacterial biosensors directly sense naphthenic acids (NA) and induce the expression of genes required to degrade or defend against NA toxicity. A panel of NA biosensors was developed where each responds to a unique subset of naphthenic acids within the complex NAFC mixtures. Bacterial gene expression is coupled to luminescence where light is produced in direct proportion to the amount of NA in a sample. A qualitative method can detect NA within 24 hrs and the quantitative method uses solid phase extraction of NA to produce an estimate of the concentration (mg/L). Biosensors confirmed orbitrap mass spectrometry monitoring of NA remediation in both mesocosm and constructed wetland water treatment systems. This simple, rapid, and specific NA detection method could be used as a high throughput NA monitoring technology during the eventual treatment and remediation of NA within oil sands tailings.



Genomics Research for Optimization of Constructed
Treatment Wetlands for Water Remediation

Connect

www.grow-genomics.ca
grow@ucalgary.ca



Click to join the Zoom meeting

